

OMNIGENE LLC

OMNIGENE LLC (USA)

SPE CARTRIDGES



HPLC COLUMN WASHING PUMP

Erida column wash pump, specially designed for LC column regeneration, one touch column cleaning function efficiently release human resource of professionals and productivity of Liquid Chromatography systems

- Double-piston pump generates stable flow and pressure to ensure the column life
- Equipped with solvent selection valve, conveniently deliver desalting and flush solvents in sequence. Also fit for other auxiliary uses
- Built-in programmable control, no need with computer or workstation. Automatically delivers two regeneration solvents in sequence
- Continuous running status detection, timely closes infusion in case of column blockage, leakage, solvent exhaustion and other incidents
- Special slim design and extreme space saving for modern laboratories
- Dedicated shelf kit, with integrated power and solvent distribution, accommodates up to 8 Erida pumps and 4 different solvents in a 30 x 30cm footprint





20 mm Aluminium GC Septa



20 mm Magnetic GC Septa



HPLC Column End Plugs




HPLC Vial



TOC Vial



PP Vial




HPLC Vial Rack



D2 Lamp



HPLC Septa



GC Graphite Ferrules



18 mm Magnetic Screw Cap



20 mm GC Crimper

Automatic Shoe Cover Machine

