

## Lab Consumable Scales & Balances Accessories



# Disposable Weighing Dishes/Boats

## Square Shape with pour-off



Disposable Weighing Dishes/Boats can be utilized for liquid or solid balance sample during weighing procedures. The smooth surface and rounded corners provide for easy removal of powdered and granular samples. Polystyrene materials flexes to allow bending of the weighing boats. Recesses at the corners facilitate controlled pour-off.

Material :PS Antistatic Color: Light Trans White

Item	Spec & PK
GTPSS041418N	Small size, Capacity: 7ml, Dim:41 x 41 x 8 mm (1.6 x 1.6 x 0.3")
GTPSSP898925N	Medium size, Capacity:100ml, Dim:89 x 89 x 25 mm (3.5 x 3.5 x 1")
GTPSSP14014025N	Large size, Capacity: 330ml, Dim:140 x 140 x 25 mm (5.5 x 5.5 x 1")

# Disposable Weighing Dishes/Boats

## Square Shape




Disposable Weighing Dishes/Boats can be utilized for liquid or solid balance sample during weighing procedures. The smooth surface and rounded corners provide for easy removal of powdered and granular samples. Polystyrene materials flexes to allow bending of the weighing boats.

Material :PS Antistatic Color: Light Trans White

Item	Spec & PK
GTPSS45458N	Small size, Capacity: 10ml, Dim: 45 x 45 x 8 mm (1.7 x 1.7 x 0.3")
GTPSS808024N	Medium size, Capacity: 100ml, Dim: 80 x 80 x 24 mm (3.2 x 3.2 x 0.9")
GTPSS13513524N	Large size, Capacity: 250ml, Dim: 135 x 135 x 24 mm (5.3 x 5.3 x 0.9")

# Disposable Weighing Dishes/Boats

## Canoe Shape



Disposable Weighing Dishes/Boats can be utilized for liquid or solid balance sample during weighing procedures. The smooth surface and rounded corners provide for easy removal of powdered and granular samples. Polystyrene materials flexes to allow bending of the weighing boats.

Material :PS Antistatic Color: Light Trans White

Item	Spec & PK
GTPSC57448N	Small size, Capacity: 20ml,Dim:57 x 44 x 8 mm (2.2 x 1.8 x 0.3")
GTPSC1378925N	Medium size, Capacity 140ml,Dim:137 x 89 x 25 mm (5.3 x 3.5 x 1")
GTPSC19112125N	Large size, Capacity: 270ml,Dim:191 x 121 x 25 mm (7.5 x 4.7 x 1")

# Disposable Weighing Dishes/Boats

## Vessel Shape



Disposable Weighing Dishes/Boats can be utilized for liquid or solid balance sample during weighing procedures. The smooth surface and rounded corners provide for easy removal of powdered and granular samples. Polystyrene materials flexes to allow bending of the weighing boats.

Material :PS Antistatic Color: Light Trans White

Item	Spec & PK
GTPSV533813N	Small size, Capacity: 15ml, Dim: 53 x 38 x 13 mm (2.1 x 1.5 x 0.5")
GTPSV1338526N	Medium size, Capacity: 120ml, Dim: 133 x 85 x 26 mm (5.2 x 3.2 x 1")
GTPSV18111026N	Large size, Capacity: 240ml, Dim: 181 x 110 x 26 mm (7.1 x 4.3 x 1")

# Aluminum Moisture Analyzer /Weighing Pans

## Moisture Analyzer Consumable

Weighing pan, is ideal for general-purpose weighing, dispensing, sampling, drying and storage, balances and moisture analyzers. Made from light gauge aluminum, no coating and is oil -free, compatible with both open and enclosed balances.



Material :Al Color: Grey

Item	Spec & PK
GTALMA1038	Capacity: 58ml, ID 103x Dep.8mm Capacity 60ml
GTALMA1285	Capacity: 90ml, Dia.128x Dep.5mm

# Aluminum Evaporating /Weighing Dishes

## Round Aluminum Dishes with Tabs



Evaporating dish, Weighing dish, is ideal for general-purpose weighing, dispensing, storage, food processing, evaporating, or dust cover. Made from light gauge aluminum, no coating and is oil free. Crimped sides for rigidity, sturdy gauge for heavy duty use. Large tab for sample numbers, flat bottom to improve sample handling, sides are angled for easy separation.

Material :Al Color: Grey

Item	Spec & PK
GTALRT5014	Small size, Capacity: 30ml, OD 50 x Dep 14mm
GTALRT6018	Medium size, Capacity: 60ml, OD 60 x Dep 18mm
GTALRT7718	Large size, Capacity: 80ml, OD 77 x Dep 18mm

# Aluminum Evaporating /Weighing Dishes

## Smooth-Walled Aluminum Dishes

Evaporating dish, Weighing dish, is ideal for general-purpose weighing, dispensing, storage, general purpose sampling and testing. Made from light gauge aluminum, no coating and is oil free. Smooth-sided, flat bottom, suit for dish contour facilitates stacking and dispensing.



Material :Al Color: Grey

Item	Spec & PK
GTALSW7218	Capacity: 60ml, Dia.72x Dep.18mm
GTALSW10319	Capacity: 135ml, Dia.103 x Dep.19mm